



Inivos®



Inivos Accessories & Consumables



Inivos®

Your partner in infection control,
delivering better care, faster

Inivos Accessories & Consumables

Save valuable time with Inivos decontamination accessories

Preparing a room for automated decontamination processes can eat into valuable staff time, and at Inivos®, we have a wide range of accessories available to speed up the preparation time, freeing up your staff to spend more time delivering care to patients.

These products have been designed for universal use with any decontamination technology, including our own Ultra-V™

and ProXcide®, and include solutions for both ultraviolet light (UV-C) and hydrogen peroxide vapour (HPV) processes.

Whether it's a system to better highlight the ongoing processes to nearby patients and staff, or keeping the room securely locked to prevent anyone walking into a space during treatment, we have a range of solutions.

ProXcin®

The ProXcide® requires our uniquely formulated solution ProXcin® to perform decontamination processes. The level of ProXcin® can be monitored via the ProXcide® and is easy replaced and disposed of once completely empty.

ProXcin® is the vapour released during the decontamination cycle which contains a low concentration of hydrogen peroxide. This low level of concentration, combined with the unique delivery method of the ProXcide®, allows for a safer and much faster re-entry time when decontaminating hospital environments and/or patient equipment.



Spill Kit

Chemical Spill Kit absorbing up to 15ltrs of chemical spill.

Kit contains:
2 x socks (60mm x 1.2m),
8 x double sided PP pads (500 x 400mm)
1 x clear plastic clip-closed carrier.



CODE	DESCRIPTION
21019	ProXcin® (4 x 2l)

CODE	DESCRIPTION
5089	Spill Kit

CapKit™

To operate a Hydrogen Peroxide Vapour system safely within your organisation, you will require a method of safely isolating the smoke detectors within spaces you are treating. The smoke detector restriction kit is a simple device for capping fire detector heads and preventing false fire alarms as a result of the vapour activating fire detection systems. A single person operation makes the CapKit™ easy to use and avoids the risk of working at height.



CODE	DESCRIPTION
5051	CapKit™ Complete
5618	Cap Cover Only
5034	TeleProp Support Pole Only

VentKit™

To operate a Hydrogen Peroxide Vapour system safely within your organisation, you will require a method of safely restricting the ventilation into the spaces you are treating. The VentKit™ is a simple device for sealing ventilation grilles by a single operator and avoiding the risks of working at height.



CODE	DESCRIPTION
5050	VentKit™ Complete 600 x 600mm
5617	Vent Cover Only 600 x 600mm
5060	VentKit™ Complete 1100 x 1100mm Inc Support Pole x 5pcs
5619	Vent Cover Only 1100 x 1100mm
5034	Vent Cover TeleProp Support Pole Only (use for 600mm & 1100mm Vent Covers)

Accessory Trolley

An easy to use and convenient way of moving the accessories you will need to complete a decontamination process, on a compact and manoeuvrable platform.

Designed to safely store:

6 x VentKits, 8 x Poles,
6 x CapKits, 1 box of ProXcin,
6 rolls of DVI tapes,
1 x Door Warning Pad



CODE	DESCRIPTION
5086	Inivos Accessory Trolley

TeleProp Pole

The TeleProp pole can be fitted to the CapKit and VentKit.

The extension range of the support pole is 1415 – 3715mm. The Teleprop comes with a quick release clip to easily extend and release after use and fitted to the base of the pole is an anti-slip rubber footing.



CODE	DESCRIPTION
5034	TeleProp Support Pole

ATI – Portasense III

Portable gas detection monitor for safe verification of air quality after decontamination using ProXcide® HPV system. Use to safely check treatment area prior to re-rentry and use.

CODE	DESCRIPTION
10509	ATI – Portasense III
10503	ATI Portasense Cell Only



DVI Sealant Tape

To operate the range of Inivos® bio-decontamination systems safely within your organisation, you will require appropriate sealant tapes to prevent leakage and visually mark the areas during the decontamination cycles for safety reasons. Our DVI Sealant Tapes have been validated for the effective containment of both Hydrogen Peroxide vapour and UV-C light when applied correctly.

CODE	DESCRIPTION
21033	DVI Sealant Tap: Low adhesion (Blue) (box of 16 rolls)
5041	DVI Sealant Tap: Medium adhesion (Red) (box of 16 rolls)
5036	DVI Sealant Tap: Maximum adhesion (Orange) (box of 8 rolls)



DoorBar™

In order to prevent access to the treatment room during a decontamination process, the DoorBar™ has been designed to swiftly secure any door. Made from force resistant extruded plastic, the DoorBar™ allows the operator to ensure the room cannot be entered, even on doors with no lock.

This is a particularly important safety consideration in the hospital environment, especially in paediatric and mental health settings.




CODE	DESCRIPTION
5055	DoorBar™

Equipment Tags

In order to operate an effective decontamination routine within your organisation, you will find it beneficial to identify equipment that has been treated with an Inivos® product or process.

These equipment tags are for single use application and manufactured from tear-proof and moisture-resistant Tyvek® material.



CODE	DESCRIPTION
5045	'Decontamination' Equipment Tags (Pack of 50)

Door Warning Signs

To operate our decontamination machines safely within your organisation, you will require a method of visually warning personnel in surrounding environments of the process and its associated hazards.

The double sided A2 Door Warning Signs are a simple and highly visible way to inform, educate and warn others of the decontamination process.



CODE	DESCRIPTION
7103	ProXcide® Door Warning Sign (Pad of 50)
5073	Ultra-V™ Door Warning Sign (Pad of 50)
20095	Ultra-V™ Connect Door Warning Sign (Pad of 50)

Door Stickers

These removable stickers can be affixed to the area around the door after a decontamination process has been completed. You will be able to tell, at a glance where and when the last successful decontamination was.

40 x 40 mm self-adhesive stickers
Supplied on a roll of 50



CODE	DESCRIPTION
5072	Ultra-V Bio Decontamination Door Sticker (Roll of 50)
7104	ProXcide Bio Decontamination Door Sticker (Roll of 50)



HERE FOR YOU AND HUNDREDS OF OTHER HOSPITALS
WORLDWIDE IN DELIVERING BETTER CARE, **FASTER**

t: +44 (0)845 270 6690 | e: customerservices@inivos.com | w: inivos.com